



Cestos de alambre de acero inoxidable

Cestos de alambre de acero inoxidable

Alambre de acero inoxidable pulido.



Code	Description	Packaging
LLG09033410	Wire baskets,st.steel,100x100x100 mm	1 pz.
LLG09033412	Wire baskets,st.steel,120x120x120 mm	1 pz.
LLG09033414	Wire baskets,st.steel,140x140x140 mm	1 pz.
LLG09033416	Wire baskets,st.steel,160x160x160 mm	1 pz.
LLG09033418	Wire baskets,st.steel,180x180x180 mm	1 pz.
LLG09033420	Wire baskets,st.steel,200x200x200 mm	1 pz.
LLG09033430	Wire baskets,st.steel,300x200x200 mm	1 pz.
LLG09033440	Wire baskets,st.steel,400x200x200 mm	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

