



## Cestos de alambre de acero inoxidable

### Cestos de alambre de acero inoxidable

Alambre de acero inoxidable pulido.



| Code        | Description                          | Packaging |
|-------------|--------------------------------------|-----------|
| LLG09033410 | Wire baskets,st.steel,100x100x100 mm | 1 pz.     |
| LLG09033412 | Wire baskets,st.steel,120x120x120 mm | 1 pz.     |
| LLG09033414 | Wire baskets,st.steel,140x140x140 mm | 1 pz.     |
| LLG09033416 | Wire baskets,st.steel,160x160x160 mm | 1 pz.     |
| LLG09033418 | Wire baskets,st.steel,180x180x180 mm | 1 pz.     |
| LLG09033420 | Wire baskets,st.steel,200x200x200 mm | 1 pz.     |
| LLG09033430 | Wire baskets,st.steel,300x200x200 mm | 1 pz.     |
| LLG09033440 | Wire baskets,st.steel,400x200x200 mm | 1 pz.     |



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5ª Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

