

**Recipientes colgantes para baños de ultrasonidos SONOREX****Bandelin****Recipientes colgantes para baños de ultrasonidos SONOREX**

Plástico sin perforar, con tapa. Para la limpieza con líquidos agresivos.



| Code        | Description                             | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG06207107 | Insert tub complete KW 10-0 for RK 510/ | 1 pz.     |
| LLG09877717 | Restoring bin KW 50 B-0from plastic wit | 1 pz.     |
| LLG09877800 | Baths, plastic KW 3                     | 1 pz.     |
| LLG09877801 | Baths, plastic, KW 5                    | 1 pz.     |
| LLG09877802 | Baths, plastic, KW 14                   | 1 pz.     |
| LLG09877803 | Baths, plastic, KW 14 B                 | 1 pz.     |
| LLG09877804 | Baths, plastic, KW 28-0                 | 1 pz.     |
| LLG09877806 | Baths, plastic, KW 50-0                 | 1 pz.     |



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5ª Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

