

**Recipientes colgantes para baños de ultrasonidos SONOREX****Bandelin****Recipientes colgantes para baños de ultrasonidos SONOREX**

Plástico sin perforar, con tapa. Para la limpieza con líquidos agresivos.



Code	Description	Packaging
LLG06207107	Insert tub complete KW 10-0 for RK 510/	1 pz.
LLG09877717	Restoring bin KW 50 B-0from plastic wit	1 pz.
LLG09877800	Baths, plastic KW 3	1 pz.
LLG09877801	Baths, plastic, KW 5	1 pz.
LLG09877802	Baths, plastic, KW 14	1 pz.
LLG09877803	Baths, plastic, KW 14 B	1 pz.
LLG09877804	Baths, plastic, KW 28-0	1 pz.
LLG09877806	Baths, plastic, KW 50-0	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5ª Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

