



Insert tubs for Sonorex ultrasonic baths

Bandelin

Insert tubs for Sonorex ultrasonic baths

Unperforated plastic, with cover. For cleaning with aggressive liquids.



Code	Description	Packaging
LLG06207107	Insert tubs for Sonorex ultrasonic baths, Type KW 10-0, For bath type RK 510/H, DT 510/H/H-RC, DL 510 H, Internal dim. (W x D x H) 242 x 182 x 135 mm	1 pz
LLG09877717	Insert tubs for Sonorex ultrasonic baths, Type KW 50 B-0, For bath type RK1050 CH, Internal dim. (W x D x H) 445 x 520 x 284 mm	1 pz
LLG09877800	Insert tubs for Sonorex ultrasonic baths, Type KW 3, For bath type RK/DT 102 H / 102 H-RC / 103 H / 100 / H, DL 102 H, Internal dim. (W x D x H) 195 x 115 x 88 mm	1 pz
LLG09877801	Insert tubs for Sonorex ultrasonic baths, Type KW 5, For bath type RK/DT 255 H/-RC, DL 255 H, Internal dim. (W x D x H) 254 x 96 x 130 mm	1 pz
LLG09877802	Insert tubs for Sonorex ultrasonic baths, Type KW 14, For bath type RK/DT 514/H, Internal dim. (W x D x H) 280 x 215 x 145 mm	1 pz
LLG09877803	Insert tubs for Sonorex ultrasonic baths, Type KW 14 B, For bath type RK/DT 514 BH/-RC, DL 514 BH, Internal dim. (W x D x H) 275 x 210 x 195 mm	1 pz
LLG09877804	Insert tubs for Sonorex ultrasonic baths, Type KW 28-0, For bath type RK 1028/H, DT 1028 H/CH, DL 1028 H, Internal dim. (W x D x H) 437 x 230 x 155 mm	1 pz
LLG09877806	Insert tubs for Sonorex ultrasonic baths, Type KW 50-0, For bath type RK 1050, Internal dim. (W x D x H) 517 x 445 x 184 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*