



Recipientes colgantes para baños de ultrasonidos SONOREX

Bandelin

Recipientes colgantes para baños de ultrasonidos SONOREX

Plástico sin perforar, con tapa. Para la limpieza con líquidos agresivos.



Code	Description	Packaging
LLG06207107	Insert tub complete KW 10-0 for RK 510/	1 pz.
LLG09877717	Restoring bin KW 50 B-0from plastic wit	1 pz.
LLG09877800	Baths, plastic KW 3	1 pz.
LLG09877801	Baths, plastic, KW 5	1 pz.
LLG09877802	Baths, plastic, KW 14	1 pz.
LLG09877803	Baths, plastic, KW 14 B	1 pz.
LLG09877804	Baths, plastic, KW 28-0	1 pz.
LLG09877806	Baths, plastic, KW 50-0	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

