



**Cestillo colgante y tapa para baños de ultrasonidos SONOREX TECHNIK**

**Bandelin**

**Cestillo colgante y tapa para baños de ultrasonidos SONOREX TECHNIK**

De acero inoxidable. Cestillos para colgar en las cubetas. Tapa protectora contra la evaporación y los ruidos.



Code	Description	Packaging
LLG09877582	Lid, stainless steel, MD 16	1 pz.
LLG09877583	Lid, stainless steel, MD 40	1 pz.
LLG09877584	Lid, stainless steel, MD 75	1 pz.
LLG09877585	Lid, stainless steel, MD 110	1 pz.
LLG09877586	Lid, stainless steel, MD 180	1 pz.
LLG09877587	Lid, stainless steel, MD 210	1 pz.
LLG09877680	Baskets, stainless steel MK 16 B for RM 1	1 pz.
LLG09877681	Baskets, stainless steel MK 40 B for RM 4	1 pz.
LLG09877682	Baskets, stainless steel MK 75 B for RM 7	1 pz.
LLG09877683	Baskets, stainless steel MK 110 for RM 11	1 pz.
LLG09877684	Baskets, stainless steel MK 180 for RM 18	1 pz.
LLG09877685	Baskets, stainless steel MK 210 for RM 21	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5ª Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

