

**Soluciones tampón técnicas****WTW****Soluciones tampón técnicas**

Identificadas mediante etiquetas de color. Tolerancia $\pm 0,03$ pH.



Code	Description	Packaging
LLG09040722	Technical buffer solutions,pH 4.01, 1 l	1 pz.
LLG09040723	Technical buffer solutions,pH 7.00, 1 l	1 pz.
LLG09040725	Technical buffer solutions,pH 10.01, 1	1 pz.
LLG09040727	Technical buffer solution 250mlph 4, TP	1 pz.
LLG09040728	Technical buffer solutions 250mlph 7,00	1 pz.
LLG09040734	Standard (DIN/NBS) buffer solution50 ml	1 pz.
LLG09040735	Standard (DIN/NBS) buffer solution50 ml	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

