

**Soluciones tampón técnicas**

Identificadas mediante etiquetas de color. Tolerancia  $\pm 0,03$  pH.



Code	Description	Packaging
LLG09040722	Technical buffer solutions,pH 4.01, 1 l	1 pz.
LLG09040723	Technical buffer solutions,pH 7.00, 1 l	1 pz.
LLG09040725	Technical buffer solutions,pH 10.01, 1	1 pz.
LLG09040727	Technical buffer solution 250mlph 4, TP	1 pz.
LLG09040728	Technical buffer solutions 250mlph 7,00	1 pz.
LLG09040734	Standard (DIN/NBS) buffer solution50 ml	1 pz.
LLG09040735	Standard (DIN/NBS) buffer solution50 ml	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

