



Soluciones tampón técnicas

WTW

Soluciones tampón técnicas

Identificadas mediante etiquetas de color. Tolerancia $\pm 0,03$ pH.



Code	Description	Packaging
LLG09040722	Technical buffer solutions,pH 4.01, 1 l	1 pz.
LLG09040723	Technical buffer solutions,pH 7.00, 1 l	1 pz.
LLG09040725	Technical buffer solutions,pH 10.01, 1	1 pz.
LLG09040727	Technical buffer solution 250mlph 4, TP	1 pz.
LLG09040728	Technical buffer solutions 250mlph 7,00	1 pz.
LLG09040734	Standard (DIN/NBS) buffer solution50 ml	1 pz.
LLG09040735	Standard (DIN/NBS) buffer solution50 ml	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

