

**Cestas de suspensión para baños de ultrasonidos SONOREX****Bandelin****Cestas de suspensión para baños de ultrasonidos SONOREX**

Acero inoxidable.

Óptima transmisión de ondas ultrasónicas. Los elementos a limpiar no deben colocarse en el fondo del tanque. K 6: Redondo, ajustable.



Code	Description	Packaging
LLG09877630	Baskets, stainless steelK 08 for RK 31	1 pz.
LLG09877631	Baskets, stainless steelK 1 for RK 52/R	1 pz.
LLG09877632	Baskets, stainless steelK 3 L for RK/Z	1 pz.
LLG09877633	Baskets, stainless steelK 3 for RK/DK 1	1 pz.
LLG09877634	Basket K 6, stainless steelfor RK 106,	1 pz.
LLG09877635	Baskets, stainless steelK 5 for 255 H	1 pz.
LLG09877636	Basket, stainless steelK6L for RK 156	1 pz.
LLG09877637	Basket, stainless steel K 6 BLfor 156 B	1 pz.
LLG09877694	Baskets, stainless steelK 14 for RK 514	1 pz.
LLG09877697	Baskets, stainless steelK 10 for RK 510	1 pz.
LLG09877698	Baskets, stainless steelK 28 for RK 1028	1 pz.
LLG09877699	Baskets, stainless steelK 50 for RK 1050	1 pz.
LLG09877700	Basket K 10 Bstainless steel,for RK 512	1 pz.
LLG09877701	Baskets, stainless steelK 14 B for 514/B	1 pz.
LLG09877702	Baskets, stainless steelK 28 C for 1028	1 pz.
LLG09877703	Baskets, stainless steelK 50C for 1050 C	1 pz.
LLG09877712	Suspension basket K 28 CAfor 6 protecti	1 pz.
LLG09877713	Restoring basket K 28 CVfor 3 full-sigh	1 pz.
LLG09877715	Restoring basket K 50 CAfor 9 brea-prot	1 pz.
LLG09877716	Restoring basket K 50 CVfor 6 full-sigh	1 pz.





**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

