

**Cestas de suspensión para baños de ultrasonidos SONOREX****Bandelin****Cestas de suspensión para baños de ultrasonidos SONOREX**

Acero inoxidable.

Óptima transmisión de ondas ultrasónicas. Los elementos a limpiar no deben colocarse en el fondo del tanque. K 6: Redondo, ajustable.



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG09877630 | Baskets, stainless steelK 08 for RK 31 | 1 pz. |
| LLG09877631 | Baskets, stainless steelK 1 for RK 52/R | 1 pz. |
| LLG09877632 | Baskets, stainless steelK 3 L for RK/Z | 1 pz. |
| LLG09877633 | Baskets, stainless steelK 3 for RK/DK 1 | 1 pz. |
| LLG09877634 | Basket K 6, stainless steelfor RK 106, | 1 pz. |
| LLG09877635 | Baskets, stainless steelK 5 for 255 H | 1 pz. |
| LLG09877636 | Basket, stainless steelK6L for RK 156 | 1 pz. |
| LLG09877637 | Basket, stainless steel K 6 BLfor 156 B | 1 pz. |
| LLG09877694 | Baskets, stainless steelK 14 for RK 514 | 1 pz. |
| LLG09877697 | Baskets, stainless steelK 10 for RK 510 | 1 pz. |
| LLG09877698 | Baskets, stainless steelK 28 for RK 1028 | 1 pz. |
| LLG09877699 | Baskets, stainless steelK 50 for RK 1050 | 1 pz. |
| LLG09877700 | Basket K 10 Bstainless steel,for RK 512 | 1 pz. |
| LLG09877701 | Baskets, stainless steelK 14 B for 514/B | 1 pz. |
| LLG09877702 | Baskets, stainless steelK 28 C for 1028 | 1 pz. |
| LLG09877703 | Baskets, stainless steelK 50C for 1050 C | 1 pz. |
| LLG09877712 | Suspension basket K 28 CAfor 6 protecti | 1 pz. |
| LLG09877713 | Restoring basket K 28 CVfor 3 full-sigh | 1 pz. |
| LLG09877715 | Restoring basket K 50 CAfor 9 brea-prot | 1 pz. |
| LLG09877716 | Restoring basket K 50 CVfor 6 full-sigh | 1 pz. |





CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

