




Tapa para baños de ultrasonidos SONOREX

Bandelin

Tapa para baños de ultrasonidos SONOREX

De acero inoxidable.

Para colocar vasos.



Code	Description	Packaging
LLG09877611	Lid,stainless steel,perforatedDE 52	1 pz.
LLG09877612	Lid,stainless steel,perforatedDE 156	1 pz.
LLG09877613	Lid,stainless steel,perforatedD 100	1 pz.
LLG09877614	Lid,stainless steel,perforatedDE 514	1 pz.
LLG09877615	Lid,stainless steel,perforatedDE 255	1 pz.
LLG09877617	Lid,stainless steel,perforatedfor 1050	1 pz.
LLG09877620	Lid,stainless steel,perforatedDE 510	1 pz.
LLG09877621	Lid,stainless steel,perforatedDE 08	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

