




**Tapa para baños de ultrasonidos SONOREX**

**Bandelin**

**Tapa para baños de ultrasonidos SONOREX**

De acero inoxidable.  
Para colocar vasos.



| Code        | Description                               | Packaging |
|-------------|---|-----------|
| LLG09877611 | Lid, stainless steel, perforated DE 52    | 1 pz.     |
| LLG09877612 | Lid, stainless steel, perforated DE 156   | 1 pz.     |
| LLG09877613 | Lid, stainless steel, perforated D 100    | 1 pz.     |
| LLG09877614 | Lid, stainless steel, perforated DE 514   | 1 pz.     |
| LLG09877615 | Lid, stainless steel, perforated DE 255   | 1 pz.     |
| LLG09877617 | Lid, stainless steel, perforated for 1050 | 1 pz.     |
| LLG09877620 | Lid, stainless steel, perforated DE 510   | 1 pz.     |
| LLG09877621 | Lid, stainless steel, perforated DE 08    | 1 pz.     |



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

