

**Soluciones electrolíticas KCL****SI Analytics****Soluciones electrolíticas KCL**

Para rellenar o sustituir el electrolito de referencia. Envasado en botellas de vidrio DURAN® de 250 ml.



Code	Description	Packaging
LLG09041033	Electrolyte solution L 4204250 ml, cali	1 pz.
LLG09041036	Electrolyte solution L 3014250 ml	1 pz.
LLG09041072	Electrolyte solutions, test solutions,2	1 pz
LLG09041092	Electrolyte solutions, test solutions,2	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

