



## Water bath RE400DB for rotary evaporators

Cole-Parmer Ltd. (Stuart)

### Water bath RE400DB for rotary evaporators

This precise digital bath has been specifically designed to be used with the Stuart rotary evaporators to give maximum efficiency for the smallest footprint. The clear LED digital display shows the actual bath temperature or the set temperature. Outer case is made from chemically resistant polypropylene.

- Digital temperature display, display resolution 1 °C, temperature control  $\pm 1$  °C
- Corrosion resistant bowl
- Easy to clean
- Integrated carry handle
- Cool touch outer case





**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

