



## Drum pumps, mains-powered

### Drum pumps, mains-powered

Lightweight and powerful pump for transferring small amounts of aqueous liquids up to 300 mPas from cans and drums with infinitely variable speed controller.

- Motor 230 V, 200 W (IP 24)
- Seal less pump tube
- Suitable accessory flow meter required
- Pump with tubing connection

Two options:

- Polypropylene version for acids and alkalis
- Stainless steel version for low-viscosity lubricants

Available as a set (9.880 310-9.880 315) with 1.5 m PVC tubing, nozzle in polypropylene, wall bracket.

**Set accessories:** Lutz nozzle, 1.5 m PVC hose 3/4", hose connection, hose clamps, wall bracket (for pump B2 Vario PP-SL Ø 32 mm).



Code	Description	Packaging
LLG09880300	Drum pumps, mains-powered, Type Pump, Immersion tube diam. 32 mm, Immersion depth 500 mm, Delivery head 7 mWS, Flow rate 75 l / min., Material PP	1 pz
LLG09880301	Drum pumps, mains-powered, Type Pump, Immersion tube diam. 32 mm, Immersion depth 700 mm, Delivery head 7 mWS, Flow rate 75 l / min., Material PP	1 pz
LLG09880302	Drum pumps, mains-powered, Type Pump, Immersion tube diam. 32 mm, Immersion depth 1000 mm, Delivery head 7 mWS, Flow rate 75 l / min., Material PP	1 pz
LLG09880303	Pump, Immersion tube diam. 28 mm, Immersion depth 500 mm, Delivery head 6.5 mWS, Flow rate 68 l / min., stainless steel	1 pz
LLG09880304	Pump, Immersion tube diam. 28 mm, Immersion depth 700 mm, Delivery head 6.5 mWS, Flow rate 68 l / min., stainless steel	1 pz
LLG09880305	Pump, Immersion tube diam. 28 mm, Immersion depth 1000 mm, Delivery head 6.5 mWS, Flow rate 68 l / min., stainless steel	1 pz
LLG09880310	Drum pumps, mains-powered, Type Pump set, Immersion tube diam. 32 mm, Immersion depth 500 mm, Delivery head 7 mWS, Flow rate 75 l / min., Material PP	1 pz
LLG09880311	Drum pumps, mains-powered, Type Pump set, Immersion tube diam. 32 mm, Immersion depth 700 mm, Delivery head 7 mWS, Flow rate 75 l / min., Material PP	1 pz
LLG09880312	Drum pumps, mains-powered, Type Pump set, Immersion tube diam. 32 mm, Immersion depth 1000 mm, Delivery head 7 mWS, Flow rate 75 l / min., Material PP	1 pz
LLG09880313	Pump set, Immersion tube diam. 28 mm, Immersion depth 500 mm, Delivery head 6.5 mWS, Flow rate 68 l / min., stainless steel	1 pz
LLG09880314	Pump set, Immersion tube diam. 28 mm, Immersion depth 700 mm, Delivery head 6.5 mWS, Flow rate 68 l / min., stainless steel	1 pz
LLG09880315	Pump set, Immersion tube diam. 28 mm, Immersion depth 1000 mm, Delivery head 6.5 mWS, Flow rate 68 l / min., stainless steel	1 pz





**CARLO ERBA Reagents S.A.S.**  
Chaussée du Vexin – 27106 Val  
de Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

