

**Bombas para barril con conexión a la red****Lutz Pumpen****Bombas para barril con conexión a la red**

Aparatos ligeros y potentes para trasegar cantidades pequeñas de medios acuosos hasta 300 mPas de bidones, toneles y recipientes de cuello estrecho con control de revoluciones.

- Motor 230V/200 W (IP 24)
- Tubo de bomba sin junta
- Se requiere un medidor de flujo accesorio adecuado
- Bomba con conexión de tubos

Dos opciones:

- Modelo de polipropileno para ácidos y bases
- Modelo de acero inoxidable para lubricantes poco viscosos

Disponible como set (9.880 310-9.880 315) con manguera de PVC de 1,5 m, boquilla de PP, soporte de pared.

**Set de accesorios:** Boquilla Lutz, manguera de PVC de 1,5 m 3/4 ", conexión de manguera, abrazaderas de manguera, soporte de pared (para bomba B2 Vario PP-SL Ø 32 mm).



Code	Description	Packaging
LLG09880300	Laboratory Vario Pump B2 PP 28-R DLwith	1 pz
LLG09880301	Laboratory Vario Pump B2 PP 28-R DLwith	1 pz
LLG09880302	Laboratory Vario Pump B2 PP 32-R SLwith	1 pz
LLG09880303	Laboratory Vario Pump B2 Niro 28-R DLwi	1 pz
LLG09880304	Laboratory Vario Pump B2 Niro 28-R DLwi	1 pz
LLG09880305	Laboratory Vario Pump B2 Niro 28-R DLwi	1 pz
LLG09880310	Pump set B2 Vario PPwith power drive, d	1 pz
LLG09880311	Pump set B2 Vario PPwith power drive, d	1 pz
LLG09880312	Pump set B2 Vario PPwith power drive, d	1 pz
LLG09880313	Pump set B2 Vario Nirowith power drive,	1 pz
LLG09880314	Pump set B2 Vario Nirowith power drive,	1 pz
LLG09880315	Pump set B2 Vario Nirowith power drive,	1 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

