



Pipetas graduadas para el cultivo de tejidos, vidrio transparente, graduación de mancha ámbar

Hirschmann

**Pipetas graduadas para el cultivo de tejidos, vidrio transparente, graduación de mancha ámbar**

Calibradas en EX(TD). Forma corta. Especiales para cultivos tisulares y para trabajar en salas limpias. Graduaciones en color ámbar resistente, marcas anulares de aforo en los puntos principales. Cero en la punta, tapón de algodón (8 mm diámetro exterior). Con pieza en la boca para el tapón de algodón.



Code	Description	Packaging
LLG06080498	Graduated pipette for tissue culture 1ml	1 pz.
LLG06080523	Graduated pipette for tissue culture 5ml	1 pz.
LLG06090340	Graduated pipette for tissue culture 2ml	1 pz.
LLG06510172	Pipette for tissue culture 10:0.1mlclae	1 pz.
LLG06510335	Pipette for tissue culture 25:0.2mlclae	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5ª Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

