




**Erlenmeyer flasks, Borosilicate glass 3.3, wide neck**

**Erlenmeyer flasks, Borosilicate glass 3.3, wide neck**

- DIN ISO 24450
- With graduations and beaded rim



Code	Description	Packaging
LLG06270181	Erlenmeyer flasks, Borosilicate glass 3.3, wide neck, 2000 ml, Neck diam. 69 mm, Bottom diam. 166 mm, Height 280 mm	10 pz
LLG09141180	Erlenmeyer flasks, Borosilicate glass 3.3, wide neck, 50 ml, Neck diam. 34 mm, Bottom diam. 51 mm, Height 85 mm	10 pz
LLG09141181	Erlenmeyer flasks, Borosilicate glass 3.3, wide neck, 100 ml, Neck diam. 34 mm, Bottom diam. 64 mm, Height 105 mm	10 pz
LLG09141182	Erlenmeyer flasks, Borosilicate glass 3.3, wide neck, 250 ml, Neck diam. 50 mm, Bottom diam. 85 mm, Height 140 mm	10 pz
LLG09141183	Erlenmeyer flasks, Borosilicate glass 3.3, wide neck, 300 ml, Neck diam. 50 mm, Bottom diam. 87 mm, Height 156 mm	10 pz
LLG09141184	Erlenmeyer flasks, Borosilicate glass 3.3, wide neck, 500 ml, Neck diam. 50 mm, Bottom diam. 105 mm, Height 175 mm	10 pz
LLG09141185	Erlenmeyer flasks, Borosilicate glass 3.3, wide neck, 1000 ml, Neck diam. 50 mm, Bottom diam. 131 mm, Height 220 mm	10 pz
LLG09141186	Erlenmeyer flasks, Borosilicate glass 3.3, wide neck, 25 ml, Neck diam. 31 mm, Bottom diam. 43 mm, Height 70 mm	10 pz
LLG09141187	Erlenmeyer flasks, Borosilicate glass 3.3, wide neck, 200 ml, Neck diam. 50 mm, Bottom diam. 79 mm, Height 131 mm	10 pz



CARLO ERBA Reagents S.A.S.  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



CARLO ERBA  
 Reagents operates with  
 a Certified Quality  
 Management System