



Tubo de PTFE

Tubo de PTFE

Politetrafluoretileno (PTFE), blanco, sin presión resistente al calor hasta 270 ° C.



Code	Description	Packaging
LLG06070190	Tubing, PTFE, Int. Ø 1mm, Ext. Ø 3mm	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

