



**Thermal printer P95DW for Gel documentation systems**

**Cleaver Scientific Ltd.**

**Thermal printer P95DW for Gel documentation systems**

The compact thermal printer with optimised heating process enables sharp printouts and clearly defined images. Thanks to reduced white space, more printouts per roll are possible.

-Resolution 325 dpi with 16-Bit data processing

-Long print format possible (100 x 450 mm)

**Scope of supply:** Thermal printer, 1 x printer paper roll





**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

