



Jeringas desechables, 3 partes, PP, estériles

Jeringas desechables, 3 partes, PP, estériles

Polipropileno con émbolo de dos partes y boquilla Luer. Sin aguja. Suministradas estériles, en bolsas selladas dentro de unas cajas.



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG06238175 | Plastipak syringes 3ml Luerlock pack200 | 200 pz |
| LLG06280595 | Plastipak® Disposable syringes PP w/120 | 120 pz |
| LLG09950300 | Plastipak syringes 20ml Luerslip pack120 | 120 pz. |
| LLG09950301 | Plastipak syringes 60ml Luerslip pack60 | 60 pz. |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

