



## Agujas desechables, PP/acero inoxidable, estériles

### Agujas desechables, PP/acero inoxidable, estériles

Con conector Luer para usar con jeringas desechables. Se suministran estériles en cajas de 100.



| Code        | Description                         | Packaging |
|-------------|-------------------------------------|-----------|
| LLG09950305 | Needles 21G x 1 1/2" green pack 100 | 100 pz.   |
| LLG09950306 | Needles 23G x 1" bleu pack 100      | 100 pz.   |
| LLG09950307 | Needles 25G x 5/8" orange pack 100  | 100 pz.   |



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5ª Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

