


**Bandeja para los baños con agitación modelos OLS/LSB****Grant****Bandeja para los baños con agitación modelos OLS/LSB**

Bandeja versátil de acero inoxidable. Diseñada para usar con una gran variedad de recipientes. Se ajusta a la capacidad máxima del recipiente gracias a que dispone de un resorte ajustable.



Code	Description	Packaging
LLG06264686	Universal tray TU26 with springs, adjus	1 pz.
LLG06264687	Universal tray TU12 with springs, adjus	1 pz.
LLG06264688	Universal tray TU18 with springs, adjus	1 pz.
LLG06264689	Test tube tray for OLS26 holds 5 x SR r	1 pz.
LLG06264690	Test tube tray for LSB12 holds 3 x SR r	1 pz.
LLG06264691	Test tube tray for LSB18 holds 5 x SR r	1 pz.
LLG06264695	Flask tray for OLS26 compatible with SC	1 pz.
LLG06264696	Flask tray for LSB12 compatible with SC	1 pz.
LLG06264697	Flask tray for LSB18 compatible with SC	1 pz.
LLG06264698	Base tray, stainless steel perforated f	1 pz.
LLG06264699	Base tray, stainless steel perforated f	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

