



Tapones para frascos de boca estrecha, series 301 / 302 / 310, LDPE

Kautex

Tapones para frascos de boca estrecha, series 301 / 302 / 310, LDPE
Colores naturales.



Code	Description	Packaging
LLG09072771	Screw cap 14 mm, LDPE for 10/20/30 ml nar	1 pz.
LLG09072773	Screw cap 18 mm, LDPE for 50/100/200 ml n	1 pz.
LLG09073314	Dropper cap with delivery tip 9 mm	1 pz.
LLG09073315	Dropper cap with delivery tip 13 mm	1 pz.
LLG09073318	Dropper cap with delivery tip 19 mm	1 pz.
LLG09073319	Cap for narrow-neck bottles LDPE with n	1 pz.
LLG09073320	Small captive cap, red PVC will be deliv	1 pz.
LLG09223110	Spray cap for 100ml bottles LDPE	1 pz.
LLG09223112	Spray cap for 250ml bottles LDPE	1 pz.
LLG09223113	Spray cap LDPE natural colored for 500ml	1 pz.
LLG09223114	Spray cap LDPE natural colored for 1000ml	1 pz.
LLG09223120	Wash heads	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

