



Tapones para frascos de boca estrecha, series 301 / 302 / 310, LDPE

Kautex

Tapones para frascos de boca estrecha, series 301 / 302 / 310, LDPE  
Colores naturales.



Code	Description	Packaging
LLG09072771	Screw cap 14 mm, LDPEfor 10/20/30 mlnar	1 pz.
LLG09072773	Screw cap 18 mm, LDPEfor 50/100/200 mln	1 pz.
LLG09073314	Dropper cap with delivery tip9 mm	1 pz.
LLG09073315	Dropper cap with delivery tip13 mm	1 pz.
LLG09073318	Dropper cap with delivery tip19 mm	1 pz.
LLG09073319	Cap for narrowl-neck bottles LDPEwith n	1 pz.
LLG09073320	Small captive cap, red PVCwill be deliv	1 pz.
LLG09223110	Spray cap for 100ml bottles LDPE	1 pz.
LLG09223112	Spray cap for 250ml bottles LDPE	1 pz.
LLG09223113	Spray cap LDPE natural coloredfor 500ml	1 pz.
LLG09223114	Spray cap LDPE natural colored for1000m	1 pz.
LLG09223120	Wash heads	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

