



Agujas para jeringas TLL, metálicas

Hamilton Central Europe SRL

Agujas para jeringas TLL, metálicas

Extraíbles, con punta Luer. Longitud de la aguja: 51 mm.



Code	Description	Packaging
LLG06089149	N NEEDLE (22/51/5) pack of 6	6 pz.
LLG06270339	N NEEDLE (22s/51/5) pack of 6	6 pz.
LLG09221657	Needles N 722s (22s/51/2) pack of 6 Atte	6 pz.
LLG09221659	Needles N 726s (26s/51/2) pack of 6 Atte	6 pz.
LLG09221672	Needles N 722 (22/51/2) pack of 6 Attent	6 pz.
LLG09221757	Needles N 722s (22s/51/3)) pack of 6 Att	6 pz.
LLG09221772	Needles N 722 pack of 6, 22/51/3	6 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

