

**Microlitre syringes, 1000 series, TLL, with metal flange**

Hamilton Central Europe SRL

**Microlitre syringes, 1000 series, TLL, with metal flange**

PTFE-coated plunger and seal, and Luer lock outlet. Gas-tight microlitre syringe without needle.



Code	Description	Packaging
LLG06083024	Microliter syringe 1025 TLL 25 ml, with	1 pz.



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

