



Microlitre syringes, 1000 series, TLL, with metal flange

Hamilton Central Europe SRL

Microlitre syringes, 1000 series, TLL, with metal flange

PTFE-coated plunger and seal, and Luer lock outlet. Gas-tight microlitre syringe without needle.



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG06083024 | Microliter syringe 1025 TLL 25 ml, with | 1 pz. |



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

