



Hidrómetros, densidad relativa (S.G.)

Hidrómetros, densidad relativa (S.G.)

Con carga granulada. Sin termómetro. Temperatura de referencia +20 °C.



Code	Description	Packaging
LLG09004006	Density hydrometers, without thermomete	1 pz.
LLG09004008	Density hydrometers,without thermometer	1 pz.
LLG09004010	Density hydrometers,without thermometer	1 pz.
LLG09004012	Density hydrometers,without thermometer	1 pz.
LLG09004014	Density hydrometers,without thermometer	1 pz.
LLG09004016	Density hydrometers,without thermometer	1 pz.
LLG09004018	Density hydrometers, without thermome/1	1 pz
LLG09004057	Density hydrometers,without thermometer	1 pz.
LLG09004060	Density hydrometers,without thermometer	1 pz.
LLG09004065	Density hydrometers,without thermometer	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

