



**Electrodos IDS para pH SenTix® 950**

Xylem Analytics Germany  
(WTW)

**Electrodos IDS para pH SenTix® 950**

**Especificaciones**

**Margen de temperatura : 0 ... 80 °C**

**Margen de pH: 0 ... 14 pH**

**Longitud/diam: 120mm/12mm**

**Electrolito:**

**Diafragma:**

**Forma de la membrana: Cilíndrica**

**Resistencia de la membrana a 25°C: 1GW**



Code	Description	Packaging
LLG07656698	pH electrode SenTix 950	1 pezzo



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

