



Electrodos IDS para pH SenTix® 950

Xylem Analytics Germany
(WTW)

Electrodos IDS para pH SenTix® 950

Especificaciones

Margen de temperatura : 0 ... 80 °C

Margen de pH: 0 ... 14 pH

Longitud/diam: 120mm/12mm

Electrolito:

Diafragma:

Forma de la membrana: Cilíndrica

Resistencia de la membrana a 25°C: 1GW



Code	Description	Packaging
LLG07656698	pH electrode SenTix 950	1 pezzo



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

