



Tubos de muestra de 20 ml para SC24/SC1/SC30, vidrio de borosilicato

Tubos de muestra de 20 ml para SC24/SC1/SC30, vidrio de borosilicato
Cristal transparente, con rosca de tornillo.



Code	Description	Packaging
LLG06252152	Sample bottles 20ml,for SC24/SC1/SC/148	148 pz
LLG07628233	Screw-lid with sealingfor sample bottle	148 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

