



Hidrómetros, densidad

Hidrómetros, densidad

DIN 12791 (y BS718). Serie M100. Can carga granulada.
Sin termómetro. Temperatura de referencia: +20 °C.



Code	Description	Packaging
LLG09004107	Density hydrometers, without thermome/1	1 pz
LLG09004108	Density hydrometers,suitable for calibr	1 pz.
LLG09004109	Density hydrometers,suitable for calibr	1 pz.
LLG09004110	Density hydrometers,suitable for calibr	1 pz.
LLG09004111	Density hydrometers,suitable for calibr	1 pz.
LLG09004112	Density hydrometers,suitable for calibr	1 pz.
LLG09004113	Density hydrometers,suitable for calibr	1 pz.
LLG09004114	Density hydrometers,suitable for calibr	1 pz.
LLG09004115	Density hydrometers,suitable for calibr	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

