



Accesorios para el analizador halógeno de humedad MB

OHAUS

Accesorios para el analizador halógeno de humedad MB



Code	Description	Packaging
LLG06260786	Needle Printer SF40A, EU with cable a/1	1 pz.
LLG06282766	Pan Holder for MB90, MB120/1	1 pz.
LLG06282767	In-Use-Cover for MB 90, MB 120/1	1 pz.
LLG06282768	Scoop, for MB90, MB120/1	1 pz.
LLG06282769	Cage, Sample for MB series/1	1 pz.
LLG06282770	Pan, Reusable, 7mm for MB series, pac/3	3 pz.
LLG06282771	Pan, Reusable, 14mm for MB series, pa/3	3 pz.
LLG09947312	Glass fibre filters 90 mm pack of 2/200	200 pezzi
LLG09947316	Anti-theft device for MB-45	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5ª Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

