


**Tubos de PFA****BOLA****Tubos de PFA**

PFA de plástico de flúor. Tienen una resistencia química casi ilimitada y se pueden usar hasta 260°C.



Code	Description	Packaging
LLG09205824	Tubing, PFA, 0.80 x 1.60 mmsupplied per	1 pz.
LLG09205825	Tubing, PFA, 1.58 x 3.18 mmsupplied per	1 pz.
LLG09205826	Tubing, PFA, 3.60 x 6.00 mmsupplied per	1 pz.
LLG09205827	Tubing, PFA, 2.00 x 4.00 mmsupplied per	1 pz.
LLG09205829	Tubing, PFA, 3.96 x 6.35 mmsupplied per	1 pz.
LLG09205830	Tubing, PFA, 4.00 x 6.00 mmsupplied per	1 pz.
LLG09205831	Tubing, PFA, 4.35 x 6.35 mmsupplied per	1 pz.
LLG09205833	Tubing, PFA, 6.00 x 8.00 mmsupplied per	1 pz.
LLG09205834	Tubing, PFA, 6.35 x 9.52 mmsupplied per	1 pz.
LLG09205835	Tubing, PFA, 6.80 x 10.00 mmsupplied pe	1 pz.
LLG09205836	Tubing, PFA, 8.00 x 10.00 mmsupplied pe	1 pz.
LLG09205839	Tubing, PFA, 10.00 x 12.00 mmsupplied p	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

