



Tapa para desecador, DURAN[®], con conexión roscada

Tapa para desecador, DURAN[®], con conexión roscada

Cristal transparente. Conexión roscada, tipo MOBILEX GL 32. DIN ISO 13130.



Code	Description	Packaging
LLG09042461	Thread lid 150mmDURAN	1 pz.
LLG09042462	Thread cover 200mm	1 pz.
LLG09042463	Thread cover 250mm	1 pz.
LLG09042467	Desiccator lid 300mmwith thread	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

