



**Tapa para desecador, DURAN® con conexión roscada**

**Tapa para desecador, DURAN® con conexión roscada**

Cristal transparente. Conexión roscada, tipo MOBILEX GL 32. DIN ISO 13130.



Code	Description	Packaging
LLG09042461	Thread lid 150mmDURAN	1 pz.
LLG09042462	Thread cover 200mm	1 pz.
LLG09042463	Thread cover 250mm	1 pz.
LLG09042467	Desiccator lid 300mmwith thread	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

