



Tapones de metal, Labocap sin mango

schuett-biotec

Tapones de metal, Labocap sin mango



| Code        | Description                              | Packaging |
|-------------|------------------------------------------|-----------|
| LLG09231334 | LABOCAP caps 12/13 mm, blackwithout/100  | 100 pz.   |
| LLG09231335 | LABOCAP caps 12/13 mm, violetwithou/100  | 100 pz.   |
| LLG09231336 | LABOCAP caps 15/16 mm, violetwithou/100  | 100 pz.   |
| LLG09231337 | LABOCAP caps 17/18 mm, blackwithout/100  | 100 pz.   |
| LLG09231338 | LABOCAP caps 17/18 mm, violetwithou/100  | 100 pz.   |
| LLG09231339 | LABOCAP caps 19/20 mm, yellowwithout/100 | 100 pz.   |
| LLG09231340 | LABOCAP caps 19/20 mm, greenwithout/100  | 100 pz.   |
| LLG09231341 | LABOCAP caps 19/20 mm, blackwithout/100  | 100 pz.   |
| LLG09231342 | LABOCAP caps 19/20 mm, violetwithou/100  | 100 pz.   |
| LLG09231369 | LABOCAP caps 9/10 mm, silverwithout/100  | 100 pz.   |
| LLG09231372 | LABOCAP caps 12/13 mm, silverwithout/100 | 100 pz.   |
| LLG09231374 | LABOCAP caps 14/15 mm, silverwithou/100  | 100 pz.   |
| LLG09231375 | LABOCAP caps 15/16 mm, silverwithou/100  | 100 pz.   |
| LLG09231377 | LABOCAP caps 17/18 mm, silverwithou/100  | 100 pz.   |
| LLG09231379 | LABOCAP caps 19/20 mm, silverwithou/100  | 100 pz.   |
| LLG09231381 | LABOCAP caps 21/23 mm, silverwithou/100  | 100 pz.   |
| LLG09231385 | LABOCAP caps 24/26 mm, silverwithou/100  | 100 pz.   |
| LLG09231388 | LABOCAP caps 28/30 mm, silverwithou/100  | 100 pz.   |
| LLG09231421 | LABOCAP caps 9/10 mm, bluewithout g/100  | 100 pz.   |
| LLG09231422 | LABOCAP caps 12/13 mm, bluewithout /100  | 100 pz.   |



| Code        | Description                             | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG09231423 | LABOCAP caps 14/15 mm, bluewithout /100 | 100 pz.   |
| LLG09231424 | LABOCAP caps 15/16 mm, bluewithout /100 | 100 pz.   |
| LLG09231425 | LABOCAP caps 17/18 mm, bluewithout /100 | 100 pz.   |
| LLG09231426 | LABOCAP caps 19/20 mm, bluewithout /100 | 100 pz.   |
| LLG09231427 | LABOCAP caps 21/23 mm, bluewithout /100 | 100 pz.   |
| LLG09231428 | LABOCAP caps 24/26 mm, bluewithout /100 | 100 pz.   |
| LLG09231429 | LABOCAP caps 28/30 mm, bluewithout /100 | 100 pz.   |
| LLG09231435 | LABOCAP-caps 12/13mm, amberw/o hand/100 | 100 pz.   |
| LLG09231441 | LABOCAP caps 9/10 mm, redwithout gr/100 | 100 pz.   |
| LLG09231442 | LABOCAP caps 12/13 mm, redwithout g/100 | 100 pz.   |
| LLG09231443 | LABOCAP caps 14/15 mm, redwithout g/100 | 100 pz.   |
| LLG09231444 | LABOCAP caps 15/16 mm, redwithout g/100 | 100 pz.   |
| LLG09231445 | LABOCAP caps 17/18 mm, redwithout g/100 | 100 pz.   |
| LLG09231446 | LABOCAP caps 19/20 mm, redwithout g/100 | 100 pz.   |
| LLG09231447 | LABOCAP caps 21/23 mm, redwithout g/100 | 100 pz.   |
| LLG09231448 | LABOCAP caps 24/26 mm, redwithout g/100 | 100 pz.   |
| LLG09231449 | LABOCAP caps 28/30 mm, redwithout g/100 | 100 pz.   |
| LLG09231460 | LABOCAP caps 12/13 mm, yellowwithou/100 | 100 pz.   |
| LLG09231461 | LABOCAP caps 15/16 mm, yellowwithou/100 | 100 pz.   |
| LLG09231462 | LABOCAP caps 17/18 mm, yellowwithou/100 | 100 pz.   |



| Code        | Description                             | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG09231463 | LABOCAP caps 12/13 mm, greenwithout/100 | 100 pz.   |
| LLG09231464 | LABOCAP caps 15/16 mm, greenwithout/100 | 100 pz.   |
| LLG09231465 | LABOCAP caps 17/18 mm, greenwithout/100 | 100 pz.   |
| LLG09231466 | LABOCAP caps 15/16 mm, blackwithout/100 | 100 pz.   |



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

