



Tapones de metal, Labocap sin mango

schuett-biotec

Tapones de metal, Labocap sin mango



Code	Description	Packaging
LLG09231334	LABOCAP caps 12/13 mm, blackwithout/100	100 pz.
LLG09231335	LABOCAP caps 12/13 mm, violetwithou/100	100 pz.
LLG09231336	LABOCAP caps 15/16 mm, violetwithou/100	100 pz.
LLG09231337	LABOCAP caps 17/18 mm, blackwithout/100	100 pz.
LLG09231338	LABOCAP caps 17/18 mm, violetwithou/100	100 pz.
LLG09231339	LABOCAP caps 19/20 mm, yellowwithou/100	100 pz.
LLG09231340	LABOCAP caps 19/20 mm, greenwithout/100	100 pz.
LLG09231341	LABOCAP caps 19/20 mm, blackwithout/100	100 pz.
LLG09231342	LABOCAP caps 19/20 mm, violetwithou/100	100 pz.
LLG09231369	LABOCAP caps 9/10 mm, silverwithout/100	100 pz.
LLG09231372	LABOCAP caps 12/13 mm, silverwithout/100	100 pz.
LLG09231374	LABOCAP caps 14/15 mm, silverwithout/100	100 pz.
LLG09231375	LABOCAP caps 15/16 mm, silverwithout/100	100 pz.
LLG09231377	LABOCAP caps 17/18 mm, silverwithout/100	100 pz.
LLG09231379	LABOCAP caps 19/20 mm, silverwithout/100	100 pz.
LLG09231381	LABOCAP caps 21/23 mm, silverwithout/100	100 pz.
LLG09231385	LABOCAP caps 24/26 mm, silverwithout/100	100 pz.
LLG09231388	LABOCAP caps 28/30 mm, silverwithout/100	100 pz.
LLG09231421	LABOCAP caps 9/10 mm, bluewithout g/100	100 pz.
LLG09231422	LABOCAP caps 12/13 mm, bluewithout /100	100 pz.



Code	Description	Packaging
LLG09231423	LABOCAP caps 14/15 mm, bluewithout /100	100 pz.
LLG09231424	LABOCAP caps 15/16 mm, bluewithout /100	100 pz.
LLG09231425	LABOCAP caps 17/18 mm, bluewithout /100	100 pz.
LLG09231426	LABOCAP caps 19/20 mm, bluewithout /100	100 pz.
LLG09231427	LABOCAP caps 21/23 mm, bluewithout /100	100 pz.
LLG09231428	LABOCAP caps 24/26 mm, bluewithout /100	100 pz.
LLG09231429	LABOCAP caps 28/30 mm, bluewithout /100	100 pz.
LLG09231435	LABOCAP-caps 12/13mm, amberw/o hand/100	100 pz.
LLG09231441	LABOCAP caps 9/10 mm, redwithout gr/100	100 pz.
LLG09231442	LABOCAP caps 12/13 mm, redwithout g/100	100 pz.
LLG09231443	LABOCAP caps 14/15 mm, redwithout g/100	100 pz.
LLG09231444	LABOCAP caps 15/16 mm, redwithout g/100	100 pz.
LLG09231445	LABOCAP caps 17/18 mm, redwithout g/100	100 pz.
LLG09231446	LABOCAP caps 19/20 mm, redwithout g/100	100 pz.
LLG09231447	LABOCAP caps 21/23 mm, redwithout g/100	100 pz.
LLG09231448	LABOCAP caps 24/26 mm, redwithout g/100	100 pz.
LLG09231449	LABOCAP caps 28/30 mm, redwithout g/100	100 pz.
LLG09231460	LABOCAP caps 12/13 mm, yellowwithou/100	100 pz.
LLG09231461	LABOCAP caps 15/16 mm, yellowwithou/100	100 pz.
LLG09231462	LABOCAP caps 17/18 mm, yellowwithou/100	100 pz.



Code	Description	Packaging
LLG09231463	LABOCAP caps 12/13 mm, greenwithout/100	100 pz.
LLG09231464	LABOCAP caps 15/16 mm, greenwithout/100	100 pz.
LLG09231465	LABOCAP caps 17/18 mm, greenwithout/100	100 pz.
LLG09231466	LABOCAP caps 15/16 mm, blackwithout/100	100 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5ª Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

