




## Hidrómetros, densidad

### Hidrómetros, densidad

Con carga granulada. Sin termómetro. Temperatura de referencia: +20 °C.



Code	Description	Packaging
LLG09004153	Density hydrometers,without thermometer	1 pz.
LLG09004154	Density hydrometers,without thermometer	1 pz.
LLG09004155	Density hydrometers,without thermometer	1 pz.
LLG09004156	Density hydrometers,without thermometer	1 pz.
LLG09004157	Density hydrometers,without thermometer	1 pz.
LLG09004158	Density hydrometers,without thermometer	1 pz.
LLG09004159	Density hydrometers,without thermometer	1 pz.
LLG09004160	Density hydrometers,without thermometer	1 pz.
LLG09004161	Density hydrometers,without thermometer	1 pz.
LLG09004162	Density hydrometers,without thermometer	1 pz.
LLG09004163	Density hydrometers,without thermometer	1 pz.
LLG09004170	Density hydrometers without thermomet/1	1 ud



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

