



Hidrómetros, densidad

Hidrómetros, densidad

Con carga granulada. Sin termómetro. Temperatura de referencia: +20 °C.



Code	Description	Packaging
LLG09004153	Density hydrometers,without thermometer	1 pz.
LLG09004154	Density hydrometers,without thermometer	1 pz.
LLG09004155	Density hydrometers,without thermometer	1 pz.
LLG09004156	Density hydrometers,without thermometer	1 pz.
LLG09004157	Density hydrometers,without thermometer	1 pz.
LLG09004158	Density hydrometers,without thermometer	1 pz.
LLG09004159	Density hydrometers,without thermometer	1 pz.
LLG09004160	Density hydrometers,without thermometer	1 pz.
LLG09004161	Density hydrometers,without thermometer	1 pz.
LLG09004162	Density hydrometers,without thermometer	1 pz.
LLG09004163	Density hydrometers,without thermometer	1 pz.
LLG09004170	Density hydrometers without thermomet/1	1 ud



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

