

**Carro de transporte para salas limpias****KEK****Carro de transporte para salas limpias**

- Material acero inoxidable (1.4301 / AISI 304)
- Superficies funcionales en KEK-Finish
- Estantes con superficie lisa o estantes con agujeros perforados ($\varnothing$ 10 mm en rejilla de 15 x 15 mm) para soportar el flujo laminar de aire
- Firmemente soldada
- Estabilizado por armadura según diseño y dimensiones
- Perfiles cuadrados cerrados en la parte superior con tapas de plástico
- 4 ruedas eléctricamente conductoras, 2 de ellas bloqueables



Code	Description	Packaging
LLG04661688	Cleanroom transport trolley with smooth	1 pz.
LLG04661689	Cleanroom transport trolley with smooth	1 pz.
LLG04661690	Cleanroom transport trolley with smooth	1 pz.
LLG04661691	Cleanroom transport trolley with perfora	1 pz.
LLG04661692	Cleanroom transport trolley with smooth	1 pz.
LLG04661693	Cleanroom transport trolley with perfora	1 pz.
LLG04661694	Cleanroom transport trolley with smooth	1 pz.
LLG04661695	Cleanroom transport trolley with perfora	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

