

**Carro de transporte para salas limpias****KEK****Carro de transporte para salas limpias**

- Material acero inoxidable (1.4301 / AISI 304)
- Superficies funcionales en KEK-Finish
- Estantes con superficie lisa o estantes con agujeros perforados (\varnothing 10 mm en rejilla de 15 x 15 mm) para soportar el flujo laminar de aire
- Firmemente soldada
- Estabilizado por armadura según diseño y dimensiones
- Perfiles cuadrados cerrados en la parte superior con tapas de plástico
- 4 ruedas eléctricamente conductoras, 2 de ellas bloqueables



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG04661688 | Cleanroom transport trolley with smooth | 1 pz. |
| LLG04661689 | Cleanroom transport trolley with smooth | 1 pz. |
| LLG04661690 | Cleanroom transport trolley with smooth | 1 pz. |
| LLG04661691 | Cleanroom transport trolley with perfora | 1 pz. |
| LLG04661692 | Cleanroom transport trolley with smooth | 1 pz. |
| LLG04661693 | Cleanroom transport trolley with perfora | 1 pz. |
| LLG04661694 | Cleanroom transport trolley with smooth | 1 pz. |
| LLG04661695 | Cleanroom transport trolley with perfora | 1 pz. |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

