



## Hidrómetros, Baumé;

### Hidrómetros, Baumé

Tipo Baumé. Con carga granulada. Sin termómetro. Temperatura de referencia: +15 °C.



Code	Description	Packaging
LLG09004305	Hydrometers,Beaumé,without thermometer0	1 pz.
LLG09004310	Hydrometers,Beaumé,without thermometer0	1 pz.
LLG09004315	Hydrometers,Beaumé,without thermometer0	1 pz.
LLG09004320	Hydrometers,Beaumé,without thermometer	1 pz.
LLG09004330	Hydrometers,Beaumé,without thermometer0	1 pz.
LLG09004340	Hydrometers,Beaumé,without thermometer0	1 pz.
LLG09004350	Hydrometers,Beaumé,without thermometer0	1 pz.
LLG09004370	Hydrometers,Beaumé,without thermometer0	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

