



## Hidrómetros, Baumé;

### Hidrómetros, Baumé

Tipo Baumé. Con carga granulada. Sin termómetro. Temperatura de referencia: +15 °C.



| Code        | Description                             | Packaging |
|-------------|---|-----------|
| LLG09004305 | Hydrometers,Beaumé,without thermometer0 | 1 pz.     |
| LLG09004310 | Hydrometers,Beaumé,without thermometer0 | 1 pz.     |
| LLG09004315 | Hydrometers,Beaumé,without thermometer0 | 1 pz.     |
| LLG09004320 | Hydrometers,Beaumé,without thermometer  | 1 pz.     |
| LLG09004330 | Hydrometers,Beaumé,without thermometer0 | 1 pz.     |
| LLG09004340 | Hydrometers,Beaumé,without thermometer0 | 1 pz.     |
| LLG09004350 | Hydrometers,Beaumé,without thermometer0 | 1 pz.     |
| LLG09004370 | Hydrometers,Beaumé,without thermometer0 | 1 pz.     |



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

