



Thermohygrometer

Thermohygrometer

Universal application thermohygrometer with outdoor probe and stand/attachment device.
Accessory wall mounting (screws or adhesive tape) are available on request.

Specifications

Measuring range: 0 ... 50 °C

switchable to °F

20 ... 99 R.H.

Resolution: 0.1 °C

1% R.H.

Dimensions (W x D x H): Approx. 65 x 2 x 97 mm



Code	Description	Packaging
LLG09235135	Thermohygrometer, Type Thermohygrometer	1 pz



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

