



Areómetro para aceites minerales

Areómetro para aceites minerales

Con termómetro de -10 a +60 °C.

Calibrable.

Temperatura de referencia +15 °C.



Code	Description	Packaging
LLG09004561	Hydrometers, mineral oil, adjustable0,6	1 pz.
LLG09004568	Hydrometers, mineral oil, adjustable0,6	1 pz.
LLG09004575	Hydrometers for mineral oils,with therm	1 pz.
LLG09004582	Hydrometers for mineral oils,with therm	1 pz.
LLG09004589	Hydrometers, mineral oil, adujstable0,	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

