



## Areómetro para aceites minerales

### Areómetro para aceites minerales

Con termómetro de -10 a +60 °C.

Calibrable.

Temperatura de referencia +15 °C.



Code	Description	Packaging
LLG09004561	Hydrometers, mineral oil, adjustable0,6	1 pz.
LLG09004568	Hydrometers, mineral oil, adjustable0,6	1 pz.
LLG09004575	Hydrometers for mineral oils,with therm	1 pz.
LLG09004582	Hydrometers for mineral oils,with therm	1 pz.
LLG09004589	Hydrometers, mineral oil, adujstable0,	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

